

Dear Healthcare Provider,

Your patient has expressed interest in the ELITONE® device. We have provided this information to support a conversation regarding this new treatment option. Additional product information is available at elitone.com

ELITONE®

The ELITONE® device is an innovative, painless, and proven treatment for female stress urinary incontinence. Worn externally, it stimulates pelvic floor muscles — to stop embarrassing leaks.

How it Works – The thin GelPad is worn externally and applied to your pelvic region. With gentle stimulation, ELITONE contracts and relaxes the pelvic floor muscles during 20-minute sessions; essentially doing the Kegel exercises for you — longer and stronger than you could on your own. **Use 4 times per week and see results in as few as 6 weeks.** It's comfortable, discreet and easy to wear. Once you experience the results, ELITONE will become an essential part of your self-care routine.

Clinically Proven Effective – In a recent clinical study, 95% of women achieved a reduction in the number of leaks per day, with an average reduction of 71% after only 6 weeks of use. Women achieved clinically significant improvements in quality of life scores and pad usage. Visit elitone.com/clinicians for additional information.

Recommended by Gynecologists – In a survey at the American College of Obstetricians and Gynecologists annual meeting, 9/10 gynecologists said they would prescribe ELITONE. The device won the “New Product of the Year” award at the 2019 My Face My Body event.

ELITONE is Ideal for:

- Busy moms
- Working women
- Postpartum care
- Aversion to vaginal devices
- Limited access to physical therapist
- Performs Kegels incorrectly

Prescribing and Insurance

ELITONE is FDA cleared and may be purchased online at elitone.com. A prescription is not required to purchase; however, a prescription may help with obtaining insurance reimbursement. The product is classified under the HCPCS code E0740 (Non-Implanted Pelvic Floor Electrical Stimulator). Coverage varies by provider. Typically, patients must have undergone and failed a 4-week documented trial of pelvic muscle exercise training (i.e. Kegels) to obtain coverage. A Detailed Written Order form is provided as an attachment. HSA and FSA accounts may be used to cover the expense. Please contact Elidah at customer@elitone.com or 978-435-4324 with questions.



Next Steps

1. Determine if ELITONE is right for you
2. [Optional] Determine if your insurer covers ELITONE
3. Visit elitone.com/order to order
4. ELITONE will be shipped directly to you
5. Begin treatment
6. Follow-up with your physician as necessary



APPLY the GelPad against your skin like a pantyliner and put your clothes back on.



ADJUST the stimulation level to provide a comfortable muscle contraction



WEAR it 20 minutes a day to get 100 pelvic floor contractions.